


Highest quality diamonds and optimal crystal orientation guarantee perfect ultrathin sections and a durable edge

Section pick-up is facilitated as the boat is horizontal allowing the water to completely fill the boat all the way round

A hydrophilic surface makes it easy to wet the cutting edge, even with a low water level





Your benefits of working with Diatome

Free customer service

Sectioning tests with biological and material research specimens of all kinds.

We send you the sections along with the surfaced sample, a report on the results obtained and a recommendation of a suitable knife.

Complete discretion when working with proprietary samples.

Re-sharpening and reworking service


A re-sharpened Diatome diamond knife demonstrates the same high quality as a new knife.

Even knives purchased in previous years can continue to be re-sharpened.

All knives can be reworked into another type of knife for no extra charge, e.g. ultra to cryo or 45° to 35°.

Exchange service

Whenever you exchange a knife we offer you a new Diatome knife at an advantageous price.



ultra 35° ultra 45°


ultra 35°

In 1989 J. C. Jésior (Ref. Jésior) demonstrated considerably reduced compression, smoother section surfaces and improved structural preservation thanks to the use of our ultra 35° knives.

In the meantime, a large number of scientists have recognized the advantages of 35° knives, in particular for sectioning Lowicryls and non-homogenous specimens, as well as non-decalcified bone, dental materials, etc.

The ultra 35° knives are perfect for sectioning relatively soft materials research specimens including metals and polymers, as well as hard specimens such as semi-conductors, superconducting oxides, catalysts, nano-crystalline ceramics, etc (Refs. Mahon, Glanvill, Swab, Quintana, Maniette, Schubert-Bischoff).

The ultra 35° knife has demonstrated its usefulness as a standard knife for the majority of applications in both biological and materials research.



The ultra 35° knife (in the triangular holder) with a cutting range of 100nm - 2µm is used for dry sectioning of epoxy or acrylic resin embedded biological samples, which need to be investigated by element analysis (Ref. Edelmann) and SIMS (Ref. Guerquin-Kern). The gliding of the sections on the dry knife surface is facilitated with the use of our Static Line II ionizer.

ultra 45°

Acknowledged as the appropriate knife angle for routine sectioning of both biological and material research specimens, it represents a balanced compromise between section quality and durability.

▶
▲
EM micrograph of an ultramicrotomed section of the anodic alumina film formed on Al-2 wt%Cu alloy. Scale bar = 100nm. Xiarong Zhou, School of Materials, University of Manchester.

◀
Ultrastructure of the roundworm *Caenorhabditis elegans*. Scale bar = 1µm. Thomas Müller-Reichert, EM Technology Development, MPI Dresden, and Kent McDonald, Electron Microscopy Laboratory, University of California, Berkeley.

◀◀
Rat muscle (Quadriceps) x 23'000
Werner Graber, Anatomisches Institut, Bern.



Specifications

ultra 35°


Cutting ranges: 30 - 200nm (standard biological applications)
10 - 30nm (Electron Spectroscopic Imaging ESI, 3D reconstruction, etc)
50nm - 1µm (alternating ultra-thin/semi-thin sectioning)
100nm - 2µm (dry sections for element analysis and SIMS)

Available sizes: 1.5mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm

ultra 45°

Cutting range: 30nm - 200nm

Available sizes: 1.5mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm



▼
Cross-section of a Si/SiGe multilayer system.
Scale bar = 10nm.
Joanne Etheridge and Selwyn Glanvill,
Electron Microscope and Microanalysis Facility,
Monash University, Melbourne.

References

- J.C. Jérior: Use of low-angle diamond knives leads to improved ultrastructural preservation of ultrathin sections. Scanning Microscopy Supplement 3, pp. 17-153, 1989. Scanning Microscopy International, Chicago (AMF O'Hare) IL 6066 USA.
- L. Edelmann: Freeze-substitution and the preservation of diffusible ions. Journal of Microscopy, Vol. 161, pp. 217-228, 1991.
- G. Mahon and T. Malis: Ultramicrotomy of Nano-crystalline Materials. Microscopy Research and Technique, Vol. 31, pp. 267-274, 1995.
- S.R. Glanvill: Ultramicrotomy of Semiconductors and Related Materials. Microscopy Research and Technique, Vol. 31, pages 267-274, 1995.
- P. Swab and R.E. Klinger: Preparation of multilayer coatings for cross-sectional Microanalysis by Ultramicrotomy. Mat. Res. Soc. Symp. Proc. Vol. 115, pages 229-234, 1989.
- P. Swab: Ultra-microscopy of Diamond Films for TEM Cross-Section Analysis. Microscopy Research and Technique, Vol. 31, pp. 308-310, 1995.
- C. Quintana: Ultramicrotomy for Cross-sections of Nanostructure. Micron Vol. 28, No. 3, pages 217-219, 1997.
- Y. Maniette: Microtomy, a convenient method for preparing TEM samples in ceramic science. Journal of Material Science Letters 9, pages 48-50, 1990.
- P. Schubert-Bischoff and T. Krist: Fast cross-sectioning technique for thin films by Ultramicrotomy. Microscopy and Microanalysis, proceedings, page 359, 1997.
- J.L. Guerquin-Kern, T.D. Wu, C. Quintana, A. Croisy: Progress in analytical imaging of the cell by dynamic secondary ion mass spectroscopy (SIMS microscopy). BBA 1724, pp. 228-238, 2005.

Diatome US
1560 Industry Road
Hatfield PA 19440
Phone 215 412 8390
Fax 215 412 8450
e-mail: sgkcck@aol.com
www.emsdiasum.com